

Prior Medical Systems B.V.

Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior OP 1010, IP 68, the multifunctional balance for quickly weighing diapers, liquids, body substances, organs.

- Capacity 3/6 or 6 /15kg
- Graduation 1gr < 3 kg > 2 , of 2 gr < 6 > 5 gr
- Plateau dimensions 220 x 150 mm, product 237(W) x 192(D) x 117 mm(H)
- Quick response, highly accurate
- Economical, 220 hr battery life
- **Class III calibrated**
- Auto off
- Source : Batteries or adaptor

For information: info@prior-medical.com