

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model PMS 3501 W, Multifunctional Baby scale to weigh baby's and infants quickly, easy and safely

- Maximum capacity 20 kg
- Graduation per **10 gram** (0 ~20 kg)
- Functions; on/off, tare (100% of capacity) and hold (if baby is moving)
- Clear display (25mm)
- Automatic on/off
- Including adaptor
- Dimensions 560 x 290 x 120 mm including 4 adjustable feet and level indicator
- WiFi data transfer
- Optional digital rod
- **Calibrated, Class III, MDD 93/42/CE amendment 2007/47/CE**

For information : info@prior-medical.com