

Prior Medical Systems B.V.  
Korfbalstraat 10  
2492 VB Den Haag

Postbus 364  
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086  
Email : [info@prior-medical.com](mailto:info@prior-medical.com)



# Prior Medical Systems



**Model Prior 5000 , Multifunctional Wheelchair and Stretcher Weighing Scale**  
for hospitals and clinics, dialysis centers and long stay care.

- Weighing capacity **300 kg**
- Graduation per **100 gram**
- Functions : BMI, Hold, Tare
- Powered by battery or adapter
- Fully electronic with 4 OIML approved loadcells
- Movable through 2 wheels at the side
- Weight : 52 kg
- Platform drive in dimensions 1200 (l) x 800 (w) mm
- Foldable handle with mounted indicator
- **OIML approved for class III, MDD 93/42/CEE amendment 2007/47/CE**

For more information email : [info@prior-medical.com](mailto:info@prior-medical.com)