

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model 6600, class III, scale for weighing the patient in the Bed

- Capacity 300 kg
- Graduation 200 grams
- Dimensions 1230 (b) x 400 (d) x 90 (h) mm
- Easily removable connection cables
- Functions: Hold and Auto Off, BMI, Tare, Save 10 tare positions
- Parts easy to connect
- Including carrying case
- Class III calibrated, MDD 93/42/CEE amendment 2007/47/CE

For information: info@prior-medical.com

