

Prior Medical Systems B.V.
Korfbalstraat 10
2492 VB Den Haag

Postbus 364
2400 AJ Alphen aan den Rijn

Tel : +31-85-0250086
Email : info@prior-medical.com



Prior Medical Systems



Prior Model BW 2200 Ultra, Multifunctional electronic personal scale with automatic BMI function and ultrasonic height meter.

- Maximum capacity 250 kg
- Classification per 50 grams < 150 kg > 100 grams
- Features; on/off, tare, hold and automatic BMI calculation
- Measuring rod 60– 210 cm, division 5 mm
- Two ultrasonic sensors, automatic calibration
- Clear display (24mm)
- Automatic on/off
- 2 Wheels for easy moving
- Adapter included
- Platform dimensions 320 x 370 mm
- including 4 adjustable feet and spirit level indication
- Class III calibrated, MDD 93/42/EC amendment 2007/47/EC
- Options : printer, voice control

For more information email : info@prior-medical.com